Voltasøjle

15.12.10

AC 3522.00

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
15.12.10

AC 3522.00

The voltaic pile is the first galvanic battery. Alessandro Volta invented the voltaic pile in the year 1800 made thereby possible the investigation of electric currents.

This model follows basically the original construction – only with modern materials. Perspex is substituted for the original wooden parts.

As the voltaic pile is a primary battery, the zinc electrodes will eventually erode – “will be eaten up”. New zinc discs are obtainable under the parts number 3522.01.

Use

The battery consists of a number of galvanic cells, each constructed like this: A zinc disc – a porous

disc moistened with an electrolyte – a copper disc. The cells are piled on top of each other, oriented the same way, i.e. the copper disc of one cell in touch with the zinc disc of the next one. A porous disc **must not** be placed **between** the cells.

As electrolyte may be used lemon juice, vinegar, salt water, copper sulfate solution or diluted sulfuric acid.

Sulfuric acid must be handled with caution! Use protective gloves.

The porous discs should not be so wet that the electrolyte is pressed out of the discs when the pile is built. It is especially unfavorable if the fluid seeps in between the cells as this will lead to increased corrosion.

When you have finished the experiments with the voltaic pile, it should be disassembled and cleaned at once. The porous discs are rinsed to remove the electrolyte before the parts are packed. It is recommended to store the zinc and copper discs separately.

Experiments

Many experiments are possible depending on the level of the pupils. The type of electrolyte as well as its concentration may be varied. The open circuit voltage can be measured directly with a digital voltmeter and by connecting the battery to different loads the internal resistance can be calculated.

When comparing measurement series you should be careful to keep those parameters constant that isn't part of the experiment: The amount of electrolyte; the force pressing the discs together; the temperature of the battery etc.

Spare parts

3522.01 Zinc discs (10 pcs.)

3522.02 Copper discs (10 pcs.)

3522.03 Porous discs (10 pcs.)

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© A/S Søren Frederiksen, Ølgod