Beskrivelse

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Description

Complete energy set containing a solar cell which supplies electricity for the production of hydrogen, which is used by the fuel cell for the production of electricity. This current can be measured / displayed in the load box.

(From left in the picture) Solar cell panel, electrolysis column, fuel cell, and load box.

Technical data

Solar cells

- Dimensions: 70 x 120 x 52 mm.
- Working temperature 10-60° C.
- Resting voltage $U_0 = 2.8$ V.
- Short-circuit current $I_{max} = 250$ mA.
- Maximum power: 0.48 W
- Voltage (at maximum power): 2.4 V.
- DC Current (max. Power): 200 mA.

PEM Electrolysis Cell

- Dimensions: 85 x 190 x 90 mm.
- Power consumption: 1.0 W.
- Operating voltage: 1.4 - 1.8 V.
- Max. current: 0.5 A.
- Hydrogen production: 3.5mL/min. (0.5A).
- Water consumption: 0.1 mL / h (at 0.3 A).


PEM Fuel cell

- Dimensions: 85 x 70 x 70 mm
- Power: 0.5 W
- Voltage: 0.4-1.0 V
- Max. current: 1.0 A.
- Hydrogen consumption: 7mL/min. (1.0A).

Load box

- Dimensions 190 x 110 x 85 mm.
- Ammeter 0-2A.
- Voltmeter 0-20 V.
- Demo motor, 0.2-3.0V.
- Demo LED, 0.6-1.55V.
- Max voltage at 1 Ω : 1.0V.
- Max voltage, other areas: 3.0V.


Demostration use

Collection and filling of distilled water

Use a spray bottle to fill demineralized or distilled water in the two lower containers. The water must cover the entire electrolysis cell when it is in operation. For normal use, top up until the water reaches the 0 mL mark.

Then mount the overflow pipes on the electrolysis tank as shown in the picture above.

Then connect the two hoses from the nozzles of the electrolysis cell to the inlet nozzles of the fuel cell (the upper ones), observe the markings H₂ and O₂ on both cells.

Finally, attach the two short hoses to the fuel cell outlet.

One end of these hoses should not be connected and left open.

In the last step, the solar panel is connected to the electrolysis cell, as shown in the picture.

Respect the polarity: Red plugs is connected with red, black with black.

If a power source other than the solar cell (to accelerate the electrolysis activity) is to be used, please note that; before connecting the cables, make sure that the current is less than 500 mA.

Turn the switch on the load box on the *Open Circuit*:

- 1) Then connect the fuel cell to the power input of the load box.
- 2) Connect the voltage input directly to the fuel cell (this is not shown in the overview picture). Turn on the load box.

**Warnings and Precautions**

Hydrogen can explode - Avoid open flames.

Wear safety goggles.

Always use distilled or demineralized water - never tap water.

Observe the minimum distance of 20 cm between the solar panel and a lamp - due to the risk of overheating.

For electrolysis with external power supply:

The voltage must never exceed 1.8 V.

The current must never exceed 500 mA.

Do not block the pipes by electrolysis to create overpressure.

Never connect the fuel cell to an external power supply.

Do not short-circuit the fuel cell (not even via the load box) - "hot spots" can occur, which damage the diaphragm.

Empty the containers of water before disposing of the electrolytic cell.