



Hygicult®

Model Chart

Hygicult® product range



Hygicult® TPC

Total Plate Count Agar containing Lecithin and Tween as neutralising agents. For monitoring total bacterial counts. Cat. No. 68010



Hygicult® E

VRBG Agar. To promote the growth of Enterobacteriaceae. Cat. No. 68012



Hygicult® E/β-GUR

VRBG Agar/β-GUR Agar. For the detection of Enterobacteriaceae and β-glucuronidase positive microbes. Cat. No. 68267



Hygicult® CF

VRB Agar. For the detection of coliforms. Cat. No. 06036



Hygicult® Y&F

Malt Agar. For the detection of fungi. Cat. No. 68013

Test procedure



Contact inoculation



Swabbing



Dipping



Drying



Incubation



Interpretation of results


Hygicult TPC



Hygicult E / Hygicult CF




Hygicult E/β-GUR (growth on β-GUR medium)



Hygicult Y&F

Yeasts



Moulds

